USER MANUAL FOR BLOOD PRESSURE MONITOR

+ HEALTH&CARE

*Blood Pressure
Mnitor*

PBPMKD550N

TABLE OF CONTENTS

IMPORTANT INFORMATION	02
CONTENTS AND DISPLAY INDICATORS	03
INTENDED USE	03
CONTRAINDICATION	03
PRODUCT DESCRIPTION	04
SPECIFICATIONS	04
NOTICE	05
SETUP AND OPERATING PROCEDURES	07
1. BATTERY LOADING	07
2. CLOCK AND DATE ADJUSTMENT	08
3. CONNECTING THE CUFF TO THE MONITOR	08
4. APPLYING THE CUFF	09
5. BODY POSTURE DURING MEASUREMENT	10
6. TAKING YOUR BLOOD PRESSURE READING	10
7. DISPLAYING STORED RESULTS	12
8. DELETING MEASUREMENTS FROM THE MEMORY	13
9. CLASSIFYING BLOOD PRESSURE AND BACKLIGHT INDICATOR	13
10. TECHNICAL ALARM DESCRIPTION	14
11. TROUBLESHOOTING (1)	15
12. TROUBLESHOOTING (2)	15
MAINTENANCE	16
EXPLANATION OF SYMBOLS	18

IMPORTANT INFORMATION

NORMAL BLOOD PRESSURE FLUCTUATION


All physical activity, excitement, stress, eating, drinking, smoking, body posture and many other activities or factors (including taking a blood pressure measurement) will influence blood pressure value. Because of this, it is mostly unusual to obtain identical multiple blood pressure readings.

Blood pressure fluctuates continually ----- day and night. The highest value usually appears in the daytime and lowest one usually at midnight. Typically, the value begins to increase at around 3:00AM, and reaches to highest level in the daytime while most people are awake and active.

Considering the above information, it is recommended that you measure your blood pressure at approximately the same time each day.

Too frequent measurements may cause injury due to blood flow interference, please always relax a minimum of 1 to 1.5 minutes between measurements to allow the blood circulation in your arm to recover. It is rare that you obtain identical blood pressure readings each time.

CONTENTS AND DISPLAY INDICATORS

INTENDED USE

Fully Automatic Electronic Sphygmomanometer is for use by medical professionals or at home and is a non-invasive blood pressure measurement system intended to measure the diastolic and systolic blood pressures and pulse rate of an adult individual by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm. The cuff circumference is limited to 22cm-48cm(approx. 8 21/32"~18 29/32").

CONTRAINDICATION

It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

PRODUCT DESCRIPTION

Based on Oscillometric methodology and silicon integrated pressure sensor, blood pressure and pulse rate can be measured automatically and non-invasively. The LCD display will show blood pressure and pulse rate. The most recent 2x30 measurements can be stored in the memory with date and time stamp. The Electronic Sphygmomanometers corresponds to the below standards: IEC 60601-1:2005/EN 60601-1:2006/AC:2010 (Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance), IEC60601-1-2:2007/EN 60601-1-2:2007 /AC:2010 (Medical electrical equipment -- Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC 80601-2-30: 2009+Cor.2010/ EN 80601-2-30: 2010(Medical electrical equipment –Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) EN 1060-1: 1995 + A1: 2002 + A2: 2009 (Non-invasive sphygmomanometers - Part 1: General requirements), EN 1060-3: 1997 + A1: 2005 + A2: 2009 (Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems), ANSI/AAMI SP-10:2002+A1:2003+A2:2006.

SPECIFICATIONS

1. Product name: Blood Pressure Monitor
2. Model: KN-550N
3. Classification: Internally powered, Type BF applied part, IPX0, No AP or APG, Continuous operation
4. Machine size: Approx. 98x98x46mm (3 27/32"x3 27/32"x1 13/16")
5. Cuff circumference: 22cm-30cm (8 21/32"-11 13/16"), 30cm-42cm (11 13/16"-16 17/32") (Optional), 42cm-48cm (16 17/32"-18 29/32") (Optional)
6. Weight: Approx. 200g (7 1/16oz) (exclude batteries and cuff)
7. Measuring method: Oscillometric method, automatic inflation and measurement
8. Memory volume: 2x30 times
9. Power source: batteries: 4x1.5V  SIZE AAA
10. Measurement range:
 - Cuff pressure: 0-300mmHg
 - Systolic: 60-260mmHg

Diastolic: 40-199mmHg
Pulse rate: 40-180 beats/minute

11. Accuracy:
 - Pressure: ± 3 mmHg
 - Pulse rate: $\pm 5\%$
12. Environmental temperature for operation: 10~40°C (50~104°F)
13. Environmental humidity for operation: $\leq 90\%$ RH
14. Environmental temperature for storage and transport: -20~55°C (-4~131°F)
15. Environmental humidity for storage and transport: $\leq 90\%$ RH
16. Environmental pressure: 80KPa-105KPa
17. Battery life: Approx 270 times
18. A list of all components belonging to the pressure measuring system, including accessories: Pump, Valve, LCD, Cuff, Sensor

These specifications are subject to change without notice.

NOTICE

1. Read all of the information in the operation guide and any other literature in the box before operating the unit.
2. Stay still, calm and rest for 5 minutes before blood pressure measurement.
3. The cuff should be placed at the same level as your heart.
4. During measurement, neither speak nor move your body and arm.
5. Measuring on same arm for each measurement.
6. Please always relax at least 1 or 1.5 minutes between measurements to allow the blood circulation in your arm to recover. Prolonged over-inflation (cuff pressure exceed 300 mmHg or maintained above 15 mmHg for longer than 3 minutes) of the bladder may cause ecchymoma of your arm.
7. Consult your physician if you have any doubt about below cases:
 - 1) The application of the cuff over a wound or inflammation diseases;
 - 2) The application of the cuff on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present;
 - 3) The application of the cuff on the arm on the side of a mastectomy;
 - 4) Simultaneously used with other monitoring medical equipments on the same limb;

- 5) Need to check the blood circulation of the user.
8. This Electronic Sphygmomanometers is designed for adults and should never be used on infants or young children. Consult your physician or other health care professionals before use on older children.
 9. Do not use this unit in a moving vehicle, This may result in erroneous measurement.
 10. Blood pressure measurements determined by this monitor are equivalent to those obtained by a trained observer using the cuff/ stethoscope auscultation method, within the limits prescribed by the American National Standard Institute, Electronic or automated sphygmomanometers.
 11. Information regarding potential electromagnetic or other interference between the blood pressure monitor and other devices together with advice regarding avoidance of such interference please see part ELECTROMAGNETIC COMPATIBILITY INFORMATION.
 12. If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a signal of will be displayed. Under this condition, the Electronic Sphygmomanometers can keep function, but the results may not be accurate, it's suggested that you consult with your physician for accurate assessment. There are 2 conditions under which the signal of IHB will be displayed:
 - 1) The coefficient of variation (CV) of pulse period >25%.
 - 2) The difference of adjacent pulse period $\geq 0.14s$, and the number of such pulse takes more than 53 percentage of the total number of pulse.
 13. Please do not use the cuff other than supplied by the manufacturer, otherwise it may bring biocompatible hazard and might result in measurement error.
 14. The monitor might not meet its performance specifications or cause safety hazard if stored or used outside the specified temperature and humidity ranges in specifications.
 15. Please do not share the cuff with other infective person to avoid cross-infection.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against

harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

16. This blood pressure monitor is verified by auscultatory method. It is recommended that you check annex B of ANSI/AAMI SP-10:2002+A1:2003+A2:2006 for details of verification method if you need.

SETUP AND OPERATING PROCEDURES

1. BATTERY LOADING

- a. Open battery cover at the back of the monitor.
- b. Load four AAA size batteries. Please pay attention to polarity.
- c. Close the battery cover. When LCD shows battery symbol *replace all batteries with new ones.*

Rechargeable batteries are not suitable for this monitor. Remove the batteries if the monitor will not be used for a month or more to avoid relevant damage of battery leakage. Avoid the battery fluid to get in your eyes. If it should get in your eyes, immediately rinse with plenty of clean water and contact a physician.


The monitor, the batteries and the cuff, must be disposed of according to local regulations at the end of their usage.

2. CLOCK AND DATE ADJUSTMENT


- a. Once you install the battery or turn off the monitor, it will enter Clock Mode, and LCD will display time and date by turns. See picture 2&2-1.

Picture 2

Picture 2-1

Picture 2-2


- b. While the monitor is in Clock Mode, pressing both the *START* and *MEMORY* button simultaneously, the month will blink at first. See picture 2-2. Press the button *START* repeatedly, the day, hour and minute will blink in turn. While the number is blinking, press the button *MEMORY* to increase the number. Keep on pressing the button *MEMORY*, the number will increase fast.
- c. You can turn off the monitor by pressing *START* button when the minute is blinking, then the time and date is confirmed.
- d. The monitor will turn off automatically after 1 minute of no operation, with the time and date unchanged.
- e. Once you change the batteries, you should readjust the time and date.

When you adjust the clock and date, the backlight is green.



3. CONNECTING THE CUFF TO THE MONITOR

Insert the cuff tubing connector into the socket in the left side of the monitor. Make certain that the connector is completely inserted to avoid air leakage during blood pressure measurements.

Avoid compression or restriction of the connection tubing during measurement, which may cause inflation error, or harmful injury due to continuous cuff pressure.

4. APPLYING THE CUFF

- Pulling the cuff end through the medial loop (the cuff is packaged like this already), turn it outward (away from your body) and tighten it and close the Velcro fastener.
- Place the cuff around a bare arm 1-2cm above the elbow joint.
- While seated, place palm upside in front of you on a flat surface such as a desk or table. Position the air tube in the middle of your arm in line with your middle finger.
- The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger between your arm and the cuff.


NOTE

- Please refer to the cuff circumference range in *SPECIFICATIONS* to make sure that the appropriate cuff is used.
- Measuring on same arm each time.
- Do not move your arm, body, or the monitor and do not move the rubber tube during measurement.


4. Stay quiet, calm for 5 minutes before blood pressure measurement.
5. Please keep the cuff clean. If the cuff becomes dirty, remove it from the monitor and clear it by hand in a mild detergent, then rinse it thoroughly in cold water. Never dry the cuff in clothes dryer or iron it. Clean the cuff after the usage of every 200 times is recommended.

5. BODY POSTURE DURING MEASUREMENT

SITTING COMFORTABLY MEASUREMENT

- a. Be seated with your feet flat on the floor, and don't cross your legs.
- b. Place palm upside in front of you on a flat surface such as a desk or table.
- c. The middle of the cuff should be at the level of the right atrium of the heart.

LYING DOWN MEASUREMENT


- a. Lie on your back.
- b. Place your arm straight along your side with your palm upside.
- c. The cuff should be placed at the same level as your heart.

6. TAKING YOUR BLOOD PRESSURE READING

- a. After applying the cuff and your body is in a comfortable position, press the "start" button. A beep is heard and all display characters are shown for self-test and the backlight is green. See picture 6. Please contact the service center if segment is missing.

Picture 6


Picture 6-1

Picture 6-2


- b. Then the current memory bank (U1 or U2) is blinking. See picture 6-1. Press *MEMORY* button to change over to other bank. See picture 6-2. Confirm your selection by pressing *START* button. The current bank can also be confirmed automatically after 5 seconds with no operation.
- c. After selecting the memory bank, the monitor starts to seek zero pressure. See picture 6-3.

Picture 6-3

Picture 6-4

Picture 6-5

- d. The monitor inflates the cuff until sufficient pressure has built up for a measurement. Then the monitor slowly releases air from the cuff and carries out the measurement. Finally the blood pressure and pulse rate will be calculated and displayed on the LCD screen separately. See picture 6-4&6-5. The irregular heartbeat symbol (if any) will be displayed on the LCD. The result will be automatically stored in the current memory bank.
- e. After measurement, the monitor will turn off automatically after 1 minute of no operation. Alternatively, you can press the *START* button to turn off the monitor manually.
- f. During measurement, you can press the *START* button to turn off the monitor manually.


Please consult a health care professional for interpretation of pressure measurements.


When the monitor displaying the measurement, the backlight shown different color according to the systolic pressure and diastolic pressure. Refer to the "CLASSIFYING BLOOD PRESSURE AND BACKLIGHT INDICATOR" section

7. DISPLAYING STORED RESULTS


a. In Clock Mode, press *MEMORY* button, the current memory bank and the amount of results in this bank will be displayed and the backlight is green. See *picture 7*. Press *START* button to change over to other bank. See *picture 7-1*. Confirm your selection by pressing *MEMORY* button. The current bank can also be confirmed automatically after 5 seconds with no operation.

Picture 7

Picture 7-1


Picture 7-2

b. After selecting the memory bank, the most recent result will be shown. See *picture 7-2*. Followed by, the blood pressure and pulse rate will be shown separately. See *picture 7-3&7-4*. The irregular heartbeat symbol (if any) will be displayed on the LCD. Press *MEMORY* button again to review the next result. See *picture 7-5*. In this way, repeatedly pressing the *MEMORY* button displays the respective results measured previously. If no result stored, LCD will show dashes as *picture 7-6*.


Rys. 7-3

Rys. 7-4

Rys. 7-5

Rys. 7-6


c. When displaying the stored results, the monitor will turn off automatically after 1 minute of no operation. You can also press the button *START* to turn off the monitor manually.

When the monitor displaying the measurement, the backlight shown different color according to the systolic pressure and diastolic pressure. Refer to the “CLASSIFYING BLOOD PRESSURE AND BACKLIGHT INDICATOR” section.

8. DELETING MEASUREMENTS FROM THE MEMORY

When any result is displaying, keeping on pressing button *MEMORY* for three seconds, all results in the memory bank will be deleted after three *BEEP*. LCD will show picture 8, and backlight is green. Press the button *MEMORY* or “start”, the monitor will turn off.

Picture 8

9. CLASSIFYING BLOOD PRESSURE AND BACKLIGHT INDICATOR

The backlight can be shown different color according to blood pressure when displaying the result. The following table figures blood pressure classification and backlight indicator. Please note that other factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration. Consult with your physician for accurate assessment, and never change your treatment by yourself.

BLOOD PRESSURE CLASSIFICATION	SBP MMHG	DBP MMHG	BACKLIGHT INDICATOR
Hypotension	≤90	or ≤60	Flashing green/amber
Optimal	91-119	and 61-79	green
Normal	120 – 129	and 80 – 84	green
High-normal	130 – 139	or 85 – 89	amber
Grade 1 Hypertension	140 – 159	or 90 – 99	flashing amber
Grade 2 Hypertension	160 – 179	or 100 – 109	red
Grade 3 Hypertension	≥ 100	or ≥ 110	flashing red

It is not intended to provide a basis of any type of rush toward emergency conditions / diagnosis based on the color scheme and that the color scheme is meant only to discriminate between the different levels of blood pressure.

10. TECHNICAL ALARM DESCRIPTION

The monitor will show *HI* or *LO* as technical alarm on LCD with no delay if the determined blood pressure (systolic or diastolic) is outside the rated range specified in part SPECIFICACIONES. In this case, you should consult a physician or check if your operation violated the instructions. The technical alarm condition (outside the rated range) is preset in the factory and cannot be adjusted or inactivated. This alarm condition is assigned as low priority according to IEC 60601-1-8. The technical alarm is non-latching and need no reset. The signal displayed on LCD will disappear automatically after about 8 seconds.

11. TROUBLESHOOTING (1)

Problem	Possible Cause	Solution
LCD Display shows abnormal result	The cuff position was not correct or it was not properly tightened	Apply the cuff correctly and try again
	Body posture was not correct during testing	Review the <i>BODY POSTURE DURING MEASUREMENT</i> sections of the instructions and re-test.
	Speaking, arm or body movement, angry, excited or nervous during testing	Re-test when calm and without speaking or moving during the test
	Irregular heartbeat (arrhythmia)	It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

12. TROUBLESHOOTING (2)

Problem	Possible Cause	Solution
LCD shows low battery symbol	Low Battery	Change the batteries

LCD shows "Er 0"	Fail to detect systolic pressure	Don't move and try again
LCD shows "Er 1"	Fail to detect systolic pressure	
LCD shows "Er 2"	Fail to detect diastolic pressure	
LCD shows "Er 3"	Pneumatic system blocked or cuff is too tight during inflation	Apply the cuff correctly and try again
LCD shows "Er 4"	Pneumatic system leakage or cuff is too loose during inflation	
LCD shows "Er 5"	Cuff pressure above 300mmHg	Measure again after five minutes. If the monitor is still abnormal, please contact the local distributor or the factory.
LCD shows "Er 6"	More than 160 seconds with cuff pressure above 15 mmHg	
LCD shows "Er 7"	EEPROM accessing error	
LCD shows "Er 8"	Device parameter checking error	
LCD shows "Er A"	Pressure sensor parameter error	
No response when you press button or load battery.	Incorrect operation or strong electromagnetic interference.	Take out batteries for five minutes, and then reinstall all batteries.

If the LCD shows low battery symbol or error message, the backlight is amber.

MAINTENANCE

Do not drop this monitor or subject it to strong impact.

Avoid high temperature and solarization. Do not immerse the monitor in water as this will result in damage to the monitor.

If this monitor is stored near freezing, allow it to acclimate to room temperature before use.

Do not attempt to disassemble this monitor.

If you do not use the monitor for a long time, please remove the batteries.

It is recommended the performance should be checked every 2 years or after repair. Please contact the service center.

Clean the monitor with a dry, soft cloth or a soft cloth squeezed well after moistened with water, diluted disinfectant alcohol, or diluted detergent.

No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairably can be supplied.

The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years, and the cuff can maintain the performance characteristics for a minimum of 1000 measurements.

It is recommended the cuff should be disinfected 2 times every week if needed (For example, in hospital or in clinique). Wipe the inner side (the side contacts skin) of the cuff by a soft cloth squeezed after moistened with Ethyl alcohol (75-90%), then dry the cuff by airing.


EXPLANATION OF SYMBOLS

PART CLASS BF come in direct contact with the patient. The cuff is part of the Class BF entering into direct contact with the patient.

ENVIRONMENT PROTECTION
– Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local Authority or retailer for recycling advice.

Temperature limitation

Consult instructions for use

Date of manufacture

European Representation

Direct current

CE0197

CE Marking

Caution

Manufacturer

SN

Product serial number

Keep product dry

Information for EMC compatibility:

The Quality of the product has been proofed and complies with the requirements of IEC 60601-1-2 (electromagnetic compliance and tests). To avoid interferences with a mobile phone we recommend to keep a distance of 3.25 m between the mobile phone and the blood pressure monitor or to shut off the mobile phone .

NOTES

Dystrybutor/Importer: **PLATINET S.A.**
PL 30-798 Kraków, ul. Christo Botewa 48
www.platinet.eu

CE0197

Andon Health Co., Ltd.

No 3 Jinping Street YaAn Road, Nankai District, Tianjin, 300190, China

Lotus Global Co., Ltd. 1 Four Seasons Terrace, West Drayton, Middlesex,
London, UB7 9GG, United Kingdom

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